

Your Practice. Your Preference. Your Productivity Solution.

## UIZION DR +

## VIZION OR + WICELESS

Viztek's ViZion DR + wireless panel provides exceptional image quality paired with advanced software features and functionality. The 14x17 panel uses enhanced Gadolinium technology and has unique features like drop notification emails.

It is ideal for use in both new systems as well as retrofit rooms. The versatile ViZion DR + wireless panel utilizes the ViZion Acquire software for image capture, manipulation, and sending. Viztek's acquire software offers useful tools to enhance the image acquisition process.

## **Panel Specifications**

Weight - 7.7lbs
Pixel Size - 150um
Plxel Array - 2304x2800
3.3lp/mm
Dimensions - 384x460x15.4mm
A/D Conversion 14 Bits

## **ViZion Software Features**



